

ELECTRO-SONIC  
**DF MOBILE**  
*Skincare Technology*



*A revolutionary  
technology,  
creating future  
beautiful-looking  
skin.*



**ENVIRON®**

a beautiful skin for a lifetime

# BASIC SKINCARE ROUTINE



STEP 1	Cleanse & Tone with your preferred Environ products.
OPTIONAL	Roll with your <b>Cosmetic</b> or <b>Gold Roll-CIT</b>
STEP 2	Apply one selected serum only <b>YOUTH+ RETINOL SERUMS</b> helps normalize the appearance for healthy-looking skin <b>YOUTH+ AVANCE ELIXIR</b> helps improve the signs of aging and sun damage <b>YOUTH+ FROWN SERUM</b> helps minimize the appearance of frown lines <b>YOUTH+ SÉRIENCE™ NIGHT SERUM</b> promotes the look of serene, luminous youthfulness <b>RADIANCE+ MELA-FADE SERUM SYSTEM</b> helps reveal an even, radiant appearance <b>MOISTURE+ HYDRATING SERUM</b> helps plump and smooth the look and feel of skin <b>COMFORT+ COLOSTRUM GEL</b> helps soothe the feel of sensitive skin
STEP 3	Dampen the provided towel and hold the side conductive panels.
STEP 4	<b>Switch the DF Mobile on</b> <b>Program One</b> - used with all serums with the exception of Colostrum Gel or a contraindication to iontophoresis. 15 second holds (indicated by beep) until target area is treated. <i>If you are contraindicated to use iontophoresis OR are using colostrum; use program three.</i> <b>Program Three</b> - 15 second holds (indicated by beep) until target area is treated.
STEP 5	Finish with balance of Environ regimen and Vitamin A step up.

