

focus on the power of
SKIN TECH+



*focus on giving your skin
the benefits of more nutrients with*
ELECTRO-SONIC DF MOBILE



ENVIRON®
a beautiful skin for a lifetime

BASIC SKINCARE ROUTINE



STEP 1	Cleanse & Tone with your preferred Environ products.
OPTIONAL	Roll with your Cosmetic or Gold Roll-CIT
STEP 2	Apply one selected serum only YOUTH+ RETINOL SERUMS helps normalize the appearance for healthy-looking skin YOUTH+ AVANCE ELIXIR helps improve the signs of aging and sun damage YOUTH+ FROWN SERUM helps minimize the appearance of frown lines YOUTH+ SERIÉENCE™ NIGHT SERUM promotes the look of serene, luminous youthfulness RADIANCE+ MELA-FADE SERUM SYSTEM helps reveal an even, radiant appearance MOISTURE+ HYDRATING SERUM helps plump and smooth the look and feel of skin COMFORT+ COLOSTRUM GEL helps soothe the feel of sensitive skin
STEP 3	Dampen the provided towel and hold the side conductive panels.
STEP 4	Switch the DF Mobile on Program One - used with all serums with the exception of Colostrum Gel or a contraindication to Iontophoresis. 15 second holds (indicated by beep) until target area is treated. <i>If you are contraindicated to use Iontophoresis OR are using colostrum; use program three.</i> Program Three - 15 second holds (indicated by beep) until target area is treated.
STEP 5	Finish with balance of Environ regimen and Vitamin A step up.